



Allows free dorsiflexion and blocks plantarflexion. Function is similar to DAFO 3. Proximal portion of DAFO 2 moves forward with the leg. The use of hinges adds complexity and cost when compared to a DAFO 3.

USE WITH PATIENTS WHO:

1. Require very tall posterior control for hyperextension of the knee.
2. Desire or need hinges for other reasons.



Hindfoot: excellent, wraparound control
Forefoot: excellent, wraparound control, hindfoot stable
Ankle: STOP plantarflexion, FREE dorsiflexion

This 2 model is shown with instep, forefoot and anterior straps, standard padding and toe rise.