
DEFENDER PEDIATRIC FOOT (FS3)

Specifically designed to endure the physical demands of a young amputee


• Ensures smooth rollover and exceptional ground compliance in a flexible design



- Patented technology and eliminated bolt holes in flex areas result in optimal strength and performance
- Virtually maintenance free design
- Rated up to 100 lbs

DEFENDER PEDIATRIC FOOT (FS3)

DEFENDER PEDIATRIC FOOT (FS3)



Clearance & Weight:

		Avg. weights	
Size	Clearance	Foot	Shell
16-18 cm	2.4 in	92g	80g
19-21 cm	3.0 in	107g	121g

Defender Stiffness Category Chart:

outegory onart.				
Weight (Ibs)	Category	Weight		
20-30	1P	9-13		
31-40	2P	14-18		
41-50	3P	19-22		
51-60	4P	23-27		
61-70	5P	28-31		
71-80	6P	32-36		
81-100	1 Adult	37-45		

REQUIRED ACCESSORIES

Foot Shell Style:

Sizes 16-18 cm- FTC-1M

Sizes 19-21 cm- FTC-2M

Choice of caucasian or ebony skin tone

OPTIONAL ACCESSORIES

- Additional Spectra Sock- Part #ACC-00-101xx (8" or 10" length)
- Additional Stiffening Bumpers-

KIT-00-119 Bumper, Defender, 40 Durometer, Sizes 16-21 cm